


STYLE 68




SMOOTH SHELL, ROUND SALINE
FILLED MAMMARY IMPLANT WITH
MODERATE PROFILE AND ANTERIOR
DIAPHRAGM VALVE. DESIGNED IN
A WIDE RANGE OF SIZES TO MEET
INDIVIDUAL PATIENT REQUIREMENTS.

JUL 28 1998

STYLE 68

SALINE-FILLED MAMMARY IMPLANT

SMOOTH SHELL, ROUND SHAPE WITH MODERATE PROFILE AND ANTERIOR DIAPHRAGM VALVE



- Position of anterior diaphragm valve is ideal for a variety of incision sites.
- Reflux valve, when attached to the fill tube, allows for easy intraoperative filling.
- Implant can be inserted through small incision sites.
- Available in twenty sizes for surgical planning convenience.

Product supplied sterile. Please see package insert for additional product information.

SIZE	CATALOG NUMBER	A DIAMETER	B PROJECTION
120-150cc	25-68121	9.0cm	3.0cm
150-180cc	25-68151	9.6cm	3.3cm
180-210cc	25-68181	10.2cm	3.4cm
210-240cc	25-68211	10.6cm	3.7cm
240-270cc	25-68241	11.1cm	3.8cm
270-300cc	25-68271	11.6cm	3.9cm
300-330cc	25-68301	11.9cm	4.1cm
330-360cc	25-68331	12.3cm	4.2cm
360-390cc	25-68361	12.7cm	4.2cm
390-420cc	25-68391	13.0cm	4.5cm
420-450cc	25-68421	13.4cm	4.5cm
450-480cc	25-68451	13.7cm	4.6cm
480-510cc	25-68481	14.1cm	4.6cm
510-540cc	25-68511	14.4cm	4.6cm
550-600cc	25-68551	14.6cm	4.9cm
600-650cc	25-68601	15.0cm	5.0cm
650-700cc	25-68651	15.2cm	5.3cm
700-750cc	25-68701	15.6cm	5.4cm
750-800cc	25-68751	15.9cm	5.6cm
800-850cc	25-68801	16.4cm	5.6cm

NOMINAL DIMENSIONS

Implant dimensions are measured at the minimum fill volume with the implant placed on a flat surface.
Note: Not every implant will conform to the dimensions given. Variations in size may occur.

FILL VOLUMES

Following the recommended fill volumes can decrease the possibility of shell wrinkling and crease-fold failure.

LIMITED WARRANTY, LIMITATION OF LIABILITY AND DISCLAIMER OF OTHER WARRANTIES

McGhan Medical warrants that reasonable care was used in the manufacture and production of this product. Because McGhan Medical has no control over the conditions of use, patient selection, surgical procedure, post-surgical stresses, or handling of the device after it leaves our possession, McGhan Medical does not warrant either a good effect or against an ill effect following its use. McGhan Medical shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from use of this product, except as provided in the McGhan Medical ConfidencePlus™ Breast Implant Replacement Program. ConfidencePlus™ provides lifetime replacement and limited financial reimbursement in the event of loss of product integrity, subject to certain conditions as fully discussed in the ConfidencePlus™ literature. McGhan Medical's sole responsibility in the event that McGhan Medical determines that the product was defective when shipped by McGhan Medical shall be replacement of the product. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law, or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for use.

CAUTION: United States Federal Law restricts this device to sale by or on the order of a licensed physician.

REGIONALLY MARKETED BY:

EUROPE
McGHAN MEDICAL B.V.
Takkebijsters 17a
4817 BL Breda, The Netherlands
Telephone: 31-765711411
Fax: 31-765712909

LATIN AMERICA
McGHAN MEDICAL MEXICO
Av Picacho Ajusco 130 DESP 701
Col Jardines en la Montaña
14210 Mexico D.F.
Telephone: 52-5-631-1155
Fax: 52-5-631-1046

ASIA-PACIFIC
McGHAN MEDICAL ASIA-PACIFIC
36a Pagoda Street
Singapore 059195
Telephone: 65-327-1570
Fax: 65-327-1571

JUL 28 1998

McGhan 
MEDICAL *Shaping Your Future...Today*

McGHAN MEDICAL CORPORATION
(800) MCGHAN1 (624-4261)
(805) 683-6761
Fax: (805) 967-5839
700 Ward Drive
Santa Barbara, CA 93111-2936 USA
www.mcghan.com

© 1998 McGhan Medical Corporation
Printed in U.S.A.

M483-D 4/98